

TangenX® SIUS® and TangenX® PRO Torque Wrench

Data Sheet

A torque wrench is required for TangenX SIUS and TangenX PRO Cassettes in the holders shown in the table below. See the TangenX Cassette User Guides for proper torque range.


- Dual Scale – calibrated dual direction at factory
- Positive lock with spring-loaded pull-down lock ring
- Fast, accurate – easy to set
- Easy-to-read laser-marked scales



Specifications

Attribute	TX019				TX026			
Torque range	30 - 250 in-lb				150 - 1000 in-lb			
Holder compatibility	TangenX SIUS TSLDI-2BMC		TangenX PRO TX002, TX003		TangenX SIUS TSPDI-4BMC		TangenX PRO TX032	
Hex socket	11/16 in				1 ¼ in			
Adapter	N/A				¾ in female to ½ in male			
Dimensions	A 11.25 in	B 1.06 in	C 0.93 in	D 1.06 in	A 16 in	B 1.43 in	C 1.20 in	D 1.40 in
Weight	0.95 lb				2.5 lb			
Accuracy	±4% CW, ±6% CCW of indicated value from 20% to 100% full scale							

Dimensions



RE_0139_07AUG2023